



Wound Care Dressings

Antimicrobial

Smith & Nephew

Acticoat® 7 Antimicrobial Dressing



- This dressing is extremely fast-acting and long-lasting, allowing it to remain on the wound for up to seven days.
- The rayon/polyester core aids in maintaining the appropriate moisture level.
- It can be cut to desired shape and size.

Size	Type	Count	Sku.No
4" x 4"	Dressing	12/Box	20101
4" x 5"	Dressing	5/Box	20141
3/4" x 12"	Rope	Each	20181
4" x 48"	Roll	Each	20501

ConvaTec

Aquacel® Ag Wound Dressing



- The first antimicrobial dressing that kills a variety of wound pathogens by providing sustained and immediate activity.
- Hydrofiber® technology provides a gelling property that when combined with wound fluid, gels on contact and locks bacteria within the gelled fibers.

Size	Count	Sku.No
2" x 2"	10/Box	3413853
4" x 4"	10/Box	3413887
4" x 4.7"	10/Box	3652336
6" x 6"	10/Box	3413895
3/4" x 18"	5/Box	3413911

Johnson & Johnson

BIOPATCH™ Antimicrobial Dressing



- This fenestrated, sterile, disk is composed of a polyurethane, absorptive foam with chlorhexidine gluconate, a well-known antiseptic agent with antimicrobial and antifungal properties.
- Ideal for use with I.V. catheters, central venous lines, arterial catheters, etc.
- Highly absorbent, stays in place and inhibits growth of bacteria for up to seven days.

Size	Count	Sku.No
3/4"	10/Box	2151
1"	10/Box	2150
1"	10/Box	3150

Antimicrobial

Healthpoint Medical

Iodosorb® Gel & Idoflex™ Pads

- Conforms to wound shape and is easily removed without disturbing delicate new tissue.
- A change in pad color from brown to yellow-gray signals a need for dressing change.
- For the treatment of chronic exuding wounds.



Size	Type	Count	Sku.No
40gm	Gel Tube	Each	2185809
4cm x 6cm	Pad	5/Box	2255958
6cm x 8cm	Pad	3/Box	2962470

Johnson & Johnson

SILVERCEL™ Antimicrobial Dressing

- Provides the broad-spectrum antimicrobial action of silver with enhanced exudate management properties of alginate technology.
- An excellent primary dressing, because the alginate forms a gel on contact with exudate.
- Do not use on third-degree burns or with patients sensitive to silver or alginates.



Size	Count	Sku.No
4 1/4" x 4 1/4"	10/Box	800404

DermaRite

SilverDerm 7™ Antimicrobial Dressing

- This dressing provides basic antimicrobial activity.
- Silver ions are released within 30 minutes of application.
- Used with chronic wounds, burns, traumatic wounds, leg ulcers, and diabetic foot ulcers.



Size	Count	Sku.No
2" x 2"	10/Box	00288
4" x 4"	10/Box	00289

Argentum Medical

Silverlon® Wound Contact Dressing

- For placement directly on wound surface, this dressing acts as a barrier and aids in reducing bacteria and fungi.
- Can be left on the wound for several days, depending on the amount of drainage.



Size	Count	Sku.No
4" x 4"	10/Box	ARWCD44

Wound Care

Dressings



Antimicrobial

Xylos Corporation

XCell® Antimicrobial Dressing



- Provides active fluid management by either absorbing or hydrating exudates, depending upon the levels.
- For use with non-infected and infected wounds; acts as a barrier and reduces the chance of infection.
- Used as a non-adherent dressing, providing less painful dressing changes.

Size	Count	Sku.No
2" x 2"	10/Box	00288
4" x 4"	10/Box	00289

Calcium Aliginate

Kendall Healthcare Co.

Curasorb™ Calcium Alginate Dressing



- Absorbs exudate from wounds and creates a protective cushioning gel.
- Once applied, it binds with exudate and helps prevent maceration.
- Provides easy and painless dressing removal, which minimizes tissue trauma and maximizes comfort.

Size	Count	Sku.No
4" x 4"	10/Box	9233
12"	5/Box	9351

DermaRite

DermaGinate® Calcium Alginate Dressing



- Readily conforms to wound bed and will not stick to healing tissue.
- May absorb up to 17 times its own weight and maintains a hydrophilic environment.
- For use with moderate to high exuding wounds.
- Creates a moist wound healing environment.

Size	Count	Sku.No
4" x 4"	10/Box	00270

Calcium Aliginate

ConvaTec

Kaltostat™ Calcium Alginate Dressing



- Soft, white, nonwoven dressing of calcium alginate fiber.
- Converts exudate to a gel/fiber mat and creates a moist, warm environment.
- May be left in place up to seven days.
- Trauma-free removal; Very absorbent.

Size	Count	Sku.No
2" x 2"	10/Box	774-224
2.8" x 4 1/2"	10/Box	68242
3" x 4 3/4"	10/Box	1600618
2gm (Rope)	5/Box	828-764

Molnlycke Health Care

Melgisorb® Alginate Dressing



- An absorbent dressing for partial to full thickness wounds with moderate to heavy exudate.
- A hydrophilic gel forms upon contact, stopping horizontal wicking and providing a moist environment suitable for wound healing.
- Appropriate for cavity and surface wounds.

Size	Count	Sku.No
4" x 4"	10/Box	251100
12 1/2" (Rope)	5/Box	253000MERC

Bertek Pharmaceuticals

Sorbsan™ Alginate Dressing



- An alginate dressing for dermal and pressure ulcers with moderate to heavy exudate
- Very absorbent and conforms easily for a perfect fit.
- Great for difficult-to-manage wounds.

Size	Count	Sku.No
2"x2"	10/Box	890-657
3"x3"	10/Box	9209
4"x4"	10/Box	1758929



Wound Care Dressings

Tubular Elastic

Western Medical

Surgitube® Tubular Gauze



- Composed of 100 percent cotton yarn knitted in a continuous, seamless tube.
- The versatile knit results in a snug fit, with a tightness depending on the tension exerted during application.
- Can be used for several layers in order to protect and secure the primary dressing.

Size	Count	SKU.No
5/8" x 5yd Size-1	Per Yard	535-187
7/8" x 5yd Size-2	Per Yard	535-195
5/8" x 50yd Size-1	Per Yard	GL-219
5/8" x 50yd Size-1	Per Yard	GL-241
7/8" x 50yd Size-2P	Per Yard	GL-242
1 1/2" x 50yd Size-3P	Per Yard	GL-243
7" x 50yd Size-3P	Per Yard	GL-225

"...established tremendous employee loyalty, helping to improve the quality of customer service."

ConvaTec

Tubifast® Retention Dressing



- Composed of a light, open weave rayon fabric that can be used as a skin covering or retention dressing for any body part.
- Holds dressing securely in place, but its light elasticity gives the patient complete freedom.
- Its unique knitted 'safety lock' does not allow this dressing to be applied too tightly.
- The open weave construction helps in avoiding maceration by allowing air to circulate over the skin.

Size	COLOR	Count	SKU.No
1 1/2" x 33'	Red	Roll	152434
2 1/8" x 33'	Green	Roll	152436
3 1/8" x 33'	Blue	Roll	152438
4 7/8" x 33'	Yellow	Roll	152440

Ointment

Healthpoint Medical

Accuzyme® Ointment

- Removes dead tissue and cleanses wounds with its enzyme and emollient combination.
- Can be used as a treatment for a variety of conditions as determined by a physician.
- Fast acting, this ointment debrides all types of necrotic tissue and is harmless to healthy tissue.



Size	Count	SKU.No
30gm	Each	2418762

Derma Sciences

Dermagran® Ointment

- This pH balanced skin protectant contains zinc-nutrient technology that provides a moist, mildly acidic environment conducive to wound healing.
- Also contains Vitamin A, Calcium, and Magnesium that aids in the healing process.
- For use with abrasions, skin tears, partial thickness pressure ulcers, and preventative skin care.



Size	Type	Count	SKU.No
4oz	Jar	Each	634-279
4oz	Tube	Each	717-728

DermaRite

DermaMed® Ointment

- This formula is pH balanced, which allows it to re-hydrate and revitalize damaged skin tissue.
- Protects and soothes chafed, burned skin, while also holding moisture in.
- Will not irritate or harm delicate tissue.



Size	Count	SKU.No
4oz	Each	00214

Wound Care

Dressings



Ointment

ConvaTec

Sensi-Care® Ointment



- This ointment seals out wetness by adhering to weepy skin, preventing further skin irritations.
- Fragile, irritated skin is soothed and protected by this special formula.
- Can be applied to clean skin in a liberal manner.

Size	Count	SkU.No
4oz	Each	325614

Surgical Trays

Busse Hospital Disposables

Busse® Dressing Change Tray

Tray contains:

- One polylined sterile field.
- One abdominal pad.
- One safety pin.
- One iodophor PVP swabstick (triples).
- Four 4" x 4" gauze sponges.
- One sharp/blunt scissors.
- One polybag with twist tie.
- One roll medical tape.
- One transfer forceps.
- One Posi-Grip™ plastic forceps



Count	SkU.No
Each	41557

Busse Hospital Disposables

Busse® E.R. Laceration Tray w/ Floor-Grade Instruments

- Contains ONE tray, CSR over-wrap, 5 inch floor-grade Webster (smooth) needle holder, 4 1/2" floor-grade straight Iris scissors, 5 inch curved floor-grade Mosquito hemostat, 1 x 2 floor-grade Adson tissue forceps, 10cc luer lock syringe, 25 gauge x 1 1/2" needle, 18 gauge x 1 1/2" needle, and polylined fenestrated drape.
- Contains TWO 2oz medicine cups and absorbent towels.
- Contains FIVE 2" x 2" gauze sponges.
- Contains TEN 4" x 4" sponges.



Count	SkU.No
20/Case	67403
20/Case	65661

Surgical Trays

Busse hospital Disposables

Busse® Sharp Debridement Tray

Tray contains:

- One pair Nitrile gloves
- One alcohol prep pad
- One PVP prep pad
- One 4" x 3" gauze sponge
- One Iris scissors
- One scalpel, #15 blade
- One metal insert forceps



Count	SkU.No
Each	20588Mercy

Busse Hospital Disposables

Busse® Suturing Tray w/ Floor-Grade Satin Finished Instruments

- Contains ONE tray, CSR over-wrap, 5"inch floor-grad Halsey needle holder with fine point satin finish, 4 1/2" floor-grade straight Iris scissors, 1 x 2 floor-grade Adson tissue forceps, polylined fenestrated drape, and absorbent towel.
- Contains TWO 2 oz medicine cups.
- Contains FIVE 2" x 2" gauze sponges.
- Contains EIGHT 4" x 4" gauze sponges.



Count	SkU.No
Each	52962

Kendal Healthcare Co.

Curity® Staple Remover Kit

- Provides an all inclusive package of quality instruments that conserves time.
- Disposable or floor grade components available.
- All kits contain one 3" x 3" gauze sponge and are sterile unless otherwise noted.



Count	SkU.No
12/Box - 48/Case	66700Mercy



Wound Care Dressings

Surgical Trays

Trinity Sterile Inc.

Trinity® Deluxe Wound Closure Tray



Tray contains:

- One 5" Baumgartner needle holder.
- One curved Kelly hemostat.
- One 5 1/4" sharp, blunt scissors.
- One 5" metal tissue forceps.
- One poly bag with twist tie.
- One polylined fenestrated drape.
- One polylined towel/drape.
- One iodophor PVP prep pad.
- One alcohol prep pad.

Count	Sku.No
20/Case	642220319

Wound Cleansers-Debriders

Healthpoint Medical

AllClenz® Wound Cleanser



- This moisturizing dermal wound cleanser washes and cleanses the wound bed, while also removing any form of debris.
- Removes debris with a wetting agent and mild surfactant.
- Will not compromise the natural healing process or harm fragile tissue.
- Creates a moist wound environment and prepares wound for dressing.

Size	Count	Sku.No
12oz	Each	020012

Bard Medical Division

Biolex® Wound Cleanser



- Helps to remove slough and debris from wound.
- Composed as a no-rinse formula with non-irritating aloe vera.
- Cleanser does not harm healthy or delicate tissue

Size	Count	Sku.No
12oz	Each	5512B

Wound Cleansers-Debriders

Carrington

Carra-Klenz™ Dermal Wound Cleanser

- This gentle, emulsifying, safe cleanser removes debris, dead tissue, and organic material from low to high exuding wounds.
- Cleanser does not interrupt the natural healing process and will not irritate the wound bed.
- Does not need to be rinsed from the wound.



Size	Count	Sku.No
6oz	Each	102060
8oz	Each	1683424
16oz	Each	102160
16oz	Each	1682855

Care-Tech Laboratories

Clinical Care® Antimicrobial Wound Cleanser

- Eschar and necrotic tissue are loosened with gentle emulsification action.
- Dressing debrides and decontaminates the wound bed and eliminates the most common organisms which inhibit wound healing.
- For use with pressure ulcers, minor burns, and severely irritated dermal areas.



Size	Count	Sku.No
12oz	24/Case	C226-12Z

Derma Sciences

Dermagran® Moisturizing Spray

- This pH balanced moisturizing spray is a skin protectant that conditions the skin before the proper dressing is applied.
- For use with Stage I pressure ulcers, abrasions, and preventative skin care.
- Composed of Vitamin A, Vitamin B6, Calcium, Magnesium, and Ionic Zinc.



Size	Count	Sku.No
4oz	Each	634-287

DermaRite

DermaKlenz® Wound Cleanser

- This non-rinsing, non-irritating, pH balanced cleanser is formulated with Zinc to promote healing.
- Easily applied and prepares the wound for the proper primary dressing.
- Cleanser aids in removing slough, organic material, debris, and blood clots without interrupting the healing process



Size	Count	Sku.No
16oz	Each	00246